



**Key Features**

- Air Cooling, Auto defrost
- Adjustable Shelves
- USB Interface

.....

**More Detail:**

## Specifications

### Technical Data

Type	Pharmacy Refrigerator
Cabinet Type	Upright
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise Level (dB)	43

### Performance

Temperature Range(°C)	2~8
-----------------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60
Power(W)	340 340
Electrical Current(A)	2.2 2.2

### Dimensions

Capacity(L/Cu.Ft)	290/10.2
Net/Gross Weight(approx)	94/116(kg) 207.2/255.7(lbs)
Interior Dimensions(W*D*H)	530*555*1080(mm) 20.9*21.9*42.5(in)

<b>Exterior Dimensions(W*D*H)</b>	665*710*1665(mm) 26.2*28.0*65.6(in)
<b>Packing Dimensions(W*D*H)</b>	735*800*1800(mm) 28.9*31.5*70.9(in)
<b>Container load (20'/40'/40'H)</b>	16/32/32

### Functions

<b>High/Low Temperature</b>	Y
<b>Remote Alarm</b>	Y
<b>Power Failure</b>	Y
<b>High/Low Voltage</b>	Y
<b>Low Battery</b>	Y
<b>Door Ajar</b>	Y

### Accessories

<b>Caster</b>	Y
<b>Foot</b>	Y
<b>Test Hole</b>	Y
<b>Shelves/Inner Doors</b>	4/-
<b>USB Interface</b>	Y

### CO2 Backup System

### LN2 Backup System

### Chart Recorder

### Others

<b>Certificate</b>	CE
--------------------	----