



PHARMACOLD

LAB-40R308

SPECIFICATIONS SHEET





TECHNICAL INFORMATION

TYPE	Vaccine Refrigerator
CABINET TYPE	Upright
COOLING TYPE	Direct Cooling
DEFROST MODE	Manual
REFRIGERANT	R290

PERFORMANCE

TEMPERATURE RANGE CELSIUS	-40°C to -20°C
ACCURACY CELSIUS	0.1°C
UNIFORMITY CELSIUS	<3.9°C

CONTROL

CONTROLLER	Microprocessor Thermostat
DISPLAY	LED

DIMENSIONS

CAPACITY (L/Cu.ft)	308
NET / GROSS WEIGHT(Approx)	225 / 251 LB
EXTERIOR DIMENSIONS(W/D/H)	27.56 / 27.17 / 75.59 in
INTERIOR DIMENSIONS(W/D/H)	21.25 / 17.71 / 50.25 in
SHELF DIMENSIONS(W/D/H)	21 / 15.5 / 6.75 in
DRAWER DIMENSIONS(W/D/H)	18.75 / 15.25 / 5.5 in

FUNCTIONS

HIGH/LOW TEMPERATURE	-
SENSOR ERROR	-

ELECTRICAL DATA

VOLTAGE (V)	115
POWER (W)	345
ELECTRICAL CURRENT (A)	3

ACCESSORIES

CASTER	Fixed
FOOT	-
PROBE HOLE	-
SHELVES/BASKET QTY	-
DOOR LOCK	Yes
REMOTE ALARM CONTACT	Yes

•**MICROCOMPUTER TEMPERATURE CONTROL**
High accuracy (0.1°C) and traceable temperature and alarm records via USB port.

•**MICRO-CHANNEL CONDENSER**
20% increased heat exchange capacity and improved compressor protection.

•**LBA FORMING TECHNOLOGY**
Better insulation, energy-saving, and larger interior storage space.

•**THREE-DIMENSIONAL SEAL**
Effective cold air retention and superior thermal insulation.

Increase the quality of your storage with Pharmacold LAB-40R308.