



# PHARMACOLD

# NAB-60

## SPECIFICATIONS SHEET



# NAB-60

## SPECIFICATIONS SHEET

### TECHNICAL INFORMATION

TYPE	Pharmacy Refrigerator
CABINET TYPE	Under Counter
COOLING TYPE	Forced Air Cooling
DEFROST MODE	Auto
REFRIGERANT	CFC-Free

### PERFORMANCE

TEMPRATURE RANGE CELSUIS	2°C to 8°C
COOLING PERFORMANCE	0.1
UNIFORMITY CELSIUS	+/-2

### CONTROL

CONTROLLER	Microprocessor Thermostat
DISPLAY	LED

### DIMENSIONS

CAPACITY (L/Cu.ft)	60 /2.1
NET / GROSS WEIGHT(Approx)	83.8/88.2 LBS
EXTERIOR DIMENSIONS(W/D/H)	20.13/22.13/26in
INTERIOR DIMENSIONS(W/D/H)	16.25/15/16.5 in

### FUNCTIONS

HIGH/LOW TEMPERATURE	Yes
POWER FAILURE	Yes
SENSOR ERROR	Yes
LOW BATTERY	Yes
DOOR AJAR	Yes

### ELECTRICAL DATA

POWER SUPPLY (V/HZ)	120/60
POWER (W)	96
ELECTRICAL CURRENT (A)	0.8

### ACCESSORIES

CASTER	No
FOOT	4
PROBE HOLE	Standard
SHELVES/BASKET QTY	1.5/-
DOOR LOCK	Standard
REMOTE ALARM CONTACT	Standard

The NAB-60 is the most versatile Pharmacold model yet. It is equipped with a reversible door which can be mounted on the left or right side. The Pharmacold NAB-60 is an asset to any medical institution or laboratory.