



PHARMACOLD

VAX-126

SPECIFICATIONS SHEET





VAX-126

SPECIFICATIONS SHEET

TECHNICAL INFORMATION

TYPE	Pharmacy Refrigerator
CABINET TYPE	Under Counter
COOLING TYPE	Forced Air Cooling
DEFROST MODE	Auto
REFRIGERANT	R600A

PERFORMANCE

TEMPRATURE RANGE CELSUIS	2°C to 8°C
ACCURACY CELSIUS	0.1°C
UNIFORMITY CELSIUS	± 1°C

CONTROL

CONTROLLER	Microprocessor Thermostat
DISPLAY	LED

DIMENSIONS

CAPACITY (L/Cu.ft)	126 / 4.4
NET / GROSS WEIGHT(Approx)	94.79 / 105.82 LBS
EXTERIOR DIMENSIONS(W/D/H)	23.5 / 24.2 /31.9 (inches)
INTERIOR DIMENSIONS(W/D/H)	19.9 / 18.4 /24.4 (inches)

FUNCTIONS

HIGH/LOW TEMPERATURE	Yes
SENSOR ERROR	Yes
POWER FAILURE	Yes
LOW BATTERY	Yes
DOOR AJAR	Yes

POWER SUPPLY (V/HZ)	115/60
POWER (W)	230
ELECTRICAL CURRENT (A)	2

ACCESSORIES

CASTER	No
FOOT	Yes
PROBE HOLE	Yes
SHELVES/DRAWER QTY	2.5
LOCK	Yes
BACKUP BATTERY	Yes

•**MICROCOMPUTER TEMPERATURE CONTROL**
Provides precise control over temperature settings.

•**TEMPERATURE ACCURACY 0.1°C**
Ensures highly accurate temperature maintenance.

Extra Features

This fridge includes a built-in USB Data Logger as well as a Remote Alarm Contact.

•**INTERNATIONAL FAMOUS COMPRESSOR**
Guarantees reliable and stable performance due to high-quality components.

•ENHANCED SAFETY

The refrigerator features a built-in battery and removable multi-layer shelves.