



# PHARMACOLD

# VAX-42

## SPECIFICATIONS SHEET



# VAX-42

## SPECIFICATIONS SHEET

### TECHNICAL INFORMATION

TYPE	Pharmacy Refrigerator
CABINET TYPE	Counter Top
COOLING TYPE	Forced Air Cooling
DEFROST MODE	Auto
REFRIGERANT	R600A

### PERFORMANCE

TEMPERATURE RANGE CELSIUS	2°C to 8°C
ACCURACY CELSIUS	0.1°C
UNIFORMITY CELSIUS	± 1°C

### CONTROL

CONTROLLER	Microprocessor Thermostat
DISPLAY	LED

### DIMENSIONS

CAPACITY (L/Cu.ft)	42 / 1.5
NET / GROSS WEIGHT(Approx)	55.11 / 59.52 LB
EXTERIOR DIMENSIONS(W/D/H)	18.9 / 17.7 /21.7 (inches)
INTERIOR DIMENSIONS(W/D/H)	15.7 / 13.6 /15.0 (inches)

### FUNCTIONS

HIGH/LOW TEMPERATURE	Yes
SENSOR ERROR	Yes
POWER FAILURE	Yes
LOW BATTERY	Yes
DOOR AJAR	Yes

POWER SUPPLY (V/HZ)	115/60
POWER (W)	115
ELECTRICAL CURRENT (A)	1

### ACCESSORIES

CASTER	No
FOOT	Yes
PROBE HOLE	Yes
SHELVES/DRAWER QTY	1/1
LOCK	Yes
BACKUP BATTERY	Yes

•**MICROCOMPUTER TEMPERATURE CONTROL**  
Provides precise control over temperature settings.

•**TEMPERATURE ACCURACY 0.1°C**  
Ensures highly accurate temperature maintenance.

•**INTERNATIONAL FAMOUS COMPRESSOR**  
Guarantees reliable and stable performance due to high-quality components.

•**ENHANCED SAFETY**

The refrigerator features a built-in battery and removable multi-layer shelves.

### Extra Features

This fridge includes a built-in USB Data Logger as well as a Remote Alarm Contact.

*Combining advanced technology and user-centric features, the Vax 42 is the ideal solution for optimal storage.*