



PHARMACOLD

VAX-756

SPECIFICATIONS SHEET



TECHNICAL INFORMATION

TYPE	Pharmacy Refrigerator
CABINET TYPE	Upright
COOLING TYPE	Forced Air Cooling
DEFROST MODE	Auto
REFRIGERANT	R600A

PERFORMANCE

TEMPERATURE RANGE CELSIUS	2°C to 8°C
ACCURACY CELSIUS	0.1°C
UNIFORMITY CELSIUS	± 1°C

CONTROL

CONTROLLER	Microprocessor Thermostat
DISPLAY	LED

DIMENSIONS

CAPACITY (L/Cu.ft)	756 / 26.7
NET / GROSS WEIGHT(Approx)	308.65 / 354.94 lbs
EXTERIOR DIMENSIONS(W/D/H)	38.2 / 31.8 /75.6 (inches)
INTERIOR DIMENSIONS(W/D/H)	34.6 / 25.0 /53.3(inches)

FUNCTIONS

HIGH/LOW TEMPERATURE	Yes
SENSOR ERROR	Yes
POWER FAILURE	Yes
LOW BATTERY	Yes
DOOR AJAR	Yes

POWER SUPPLY (V/HZ)	115/60
POWER (W)	320
ELECTRICAL CURRENT (A)	1.8

ACCESSORIES

CASTER	Yes
FOOT	Yes
PROBE HOLE	Yes
SHELVES/DRAWER QTY	10/0
LOCK	Yes
BACKUP BATTERY	Yes

•**MICROCOMPUTER TEMPERATURE CONTROL**
Provides precise control over temperature settings.

•**TEMPERATURE ACCURACY 0.1°C**
Ensures highly accurate temperature maintenance.

•**INTERNATIONAL FAMOUS COMPRESSOR**
Guarantees reliable and stable performance due to high-quality components.

•**ENHANCED SAFETY**

The refrigerator features a built-in battery and removable multi-layer shelves.

Extra Features

This fridge includes a back up USB Data Logger as well as a Remote Alarm Contact.